


iCam Can Be Used By Anyone, Including Patients!

Why iCam?

- Brilliant high-quality images!
- Simple to install and easy to use
- Works with all major software
- 1-Year Warranty with first-class installation service and technical support



Co-Inventor of the intraoral camera

In 1989 we partnered with Panasonic and ETS Groux Optical to develop a revolutionary handpiece for the dental industry; This design was the pre-cursor to all modern intraoral handpieces.

Seeing is Believing!

This Direct USB camera is ideal to motivate and educate your patients. Its image quality rivals that of more expensive cameras and the Macro lens cap boosts the magnification to 50x for superb details.



Randy Ahlquist
22424 E Weaver Drive
Aurora, CO 80016
(303) 888-5500
www.integritydental.biz
randy@integritydental.biz

  Windows 10, 8,7, Vista, XP

DENTERPRISE
INTERNATIONAL, INC

The camera is easy to use and works with most imaging software. Unlike most inexpensive cameras which have a short life span, render a mediocre image, and whose capture button does not save an image in your software, iCam performs admirably, delivering an excellent image. With its programmable capture button it is compatible with all major dental imaging software.

This quality camera comes with a 1-year warranty. First-class installation service and technical support are available. iCam is worth a closer look! Call Us Today!


Technical Specifications

- Angle of Lens: 90 Degrees
- Camera Focus: Fixed
- Depth of Field: 40mm
- Lighting Source: 6 Light Emitting Diodes
- Warranty: 1 Year Parts and Labor
- Footprint: Less than 6 Sq inch
- Interface: Universal Serial Bus (USB)
- Resolution: 1.3MP



Compatible Software

DEXIS10 ...AND MANY MORE



Randy Ahlquist
22424 E Weaver Drive
Aurora, CO 80016
(303) 888-5500
www.integritydental.biz
randy@integritydental.biz

Windows 10, 8,7, Vista, XP

