



Compare! <sup>Sug.</sup> Retail

\$5,500

-HELD

## Neck Strap Included!





**AMBIDEXTROUS** 

Indicator switch flips the active trigger from right to left hand firing.



Large bright easy to read display.

## **High-Capacity Battery**



New Fast charge internal battery - 300-500 shots per battery charge (dependent on exposure time)!



Keep your office running even if you forget to put your MaxRay on the charger. Or if the internal battery fails!





MaxRay

- Includes: Unit with Internal Battery, Back-up Battery, Charger, Charging Pad, Neck Strap, Case and Instructions Black (BX-3000K) White (BX-3000W)

MaxRay	Spec.
КVр	65
mA	1.7
Focal Spot	0.4mm
Filtration	2.3mmAl
Volts	12-18.0
Style	Camera





Compare! Sug. Retail

## MaxRay COCOON • Very light - under 5 pounds

- Powerful 70 Kvp •
- High capacity battery 500 shots
- Easy to use innovative display and voice prompts



Innovative display - full function controls. With voice prompting - yes the Cocoon TALKS!

Easy 1 handed operation. Unit easily mounts to charging dock with one hand.





MaxRay Cocoon - Includes: Unit with Internal Battery, Back-up Battery, Charging Cradle, Neck Strap, Power Adapter, Case and Instructions (BX-7017)

MaxRay Cocoon	Spec.
КVр	70
mA	1.7
Focal Spot	0.4mm
Filtration	2.3mmAl
Volts	
Style	Gun Type

Vector R&D, Inc. | 6824 19th St #230 | University Place, WA 98466 | 877.883.7455 | IridiumDental.com