



American-Made Quality Meets State-of-the-Art Design

American Made

All chairs, lights, delivery systems and supporting equipment are engineered and manufactured under our own roof in Hillsboro, Oregon.

Customizable

Our products are designed to complement any office design, space constrictions, functional requirements and price point.

World-Class Service

We give you access to a team of technicians and the manufacturing team itself for fast, efficient service and peace of mind.

EARNs ITS LIVING DOUBLES



AS A DENTAL CHAIR. AS A PATIENT PLEASER.

engle
dental systems





For Over 50 Years,

Engle has built its reputation on quality, ingenuity, and integrity. With industry-leading engineering and customer service that dentists rely on, Engle products are American-made, designed and manufactured in Hillsboro, Oregon.

Engle Dental is the only company that makes the fully programmable, 10.25" forward/backward hydraulic Traverse Chair™. The traverse feature allows dentists to position themselves behind the patient comfortably without sacrificing space. Since the seat moves forward as the back of the chair reclines, dentists can stay put and let the patient come to them. The Traverse Chair™ keeps the patient's head/oral cavity in the same vertical plane in the operatory (front to back) regardless of the chair position.

Engineered for ergonomic efficiency, the Engle Traverse Chair™ allows for unparalleled dental team access and comfort. With an anatomically correct back pivot point and synchronized back and toe board recline, the Traverse Chair™ has undergone extensive comfort testing to provide maximum comfort in procedures lasting 45 minutes or more.

With options like heat and massage, lumbar support, a double articulating headrest, and Ultraleather® upholstery, dentists can focus on patient care knowing that Engle has created a level of comfort that helps patients relax.

Table of Contents

DELIVERY SYSTEMS.....	10
- Mobile Delivery	10
- 12 O’Clock Delivery System	15
- Wall Mounted Delivery	18
- Hygiene Delivery	20
- Assistant Packages	21
PACKAGES.....	22
• 360 Packages.....	22
- Traverse	
- Standard	
• 320 Packages.....	24
- Traverse	
- Standard	
• 310 Packages.....	26
- Traverse	
- Standard	
LIGHTS.....	28
- LED Operatory Lights	28
- Monitor Assembly	29
CHAIRS.....	30
UPHOLSTERY.....	32
STOOLS.....	34
TECHNICAL DIMENSIONS.....	36



Contact Info

TECHNICAL SUPPORT

TOLL FREE: 1-800-547-1906

support@engledental.com

Available: 6:00am-3:30pm PST Monday-Friday

CUSTOMER SERVICE

TOLL FREE: 1-800-547-1906

info@engledental.com

SHIPPING ADDRESS

7205 NE Evergreen Pkwy, Suite 100
Hillsboro, OR 97124

SHIPPING INFORMATION

The authorized dealer is responsible for providing a correct shipping address and details for delivery. Products are shipped F.O.B. Hillsboro, Oregon, USA. All items should be inspected upon receipt. Any shipping damage must be filed by the receiver with the carrier.

RETURN POLICY

No products will be accepted for return without prior manufacturer authorization via RMA number. RMA can be generated from Engle's technical support team. RMA number must be marked on the outside of boxes being returned, otherwise they will be refused and returned. All returned goods are subject to a 20% restocking fee. Returned goods must be received within 30 days of RMA number. All shipping charges for return of products are the sole responsibility of the dealer.

CANCELLATION POLICY

Changes or cancellations must be received by Engle Dental within 72 hours after order acknowledgment is received. Any changes or cancellation beyond that time period will be subject to additional charges or delays to order. Contact Engle Customer Service for additional questions to specific orders. For products with custom colors, including but not limited to, dental chair upholstery, stools, dental furniture, and table tops, no return will be accepted and no credits will be available. Equipment and Parts are sold exclusively through authorized Engle Dealers. Only authorized Engle dealers are permitted to sell, service Engle Equipment and Parts. Contact Engle Customer Service for a list of Authorized Equipment Dealers.

WARRANTY

All Products except stools, sold to and installed by authorized Engle Dental Systems dealers are guaranteed against defects in materials and workmanship for a period of five (5) years after date of purchase. Stools are guaranteed against defects in materials and workmanship for a period of two (2) years after date of purchase. Upholstery is guaranteed for one (1) year after date of purchase. Handpiece illumination systems are guaranteed for six (6) months after date of purchase. This warranty does not cover light bulbs of any kind. Faulty parts as determined by dealer service representative, will be repaired or replaced; no claim for labor or consequential damages will be allowed. Invoices for replacement parts may be subject to adjustment after the manufacturer has inspected the alleged defective parts. The warranty does not cover normal wear, stains, cuts, or scratches of upholstery surface finishes, or parts sold to O.E.M. customers. This warranty does not cover damage resulting from the use of cleaning, disinfecting, or sterilization chemicals and processes. The warranty does not cover damage from improper installation or maintenance, accident, misuse, or accidents in shipping and handling. In the event of shipping damage, the purchaser must initiate all claims against the freight carrier in writing within twenty-four (24) hours after the damaged items are received. The claim is the responsibility of the customer.



Ordering Information

ORDERING INFORMATION

When ordering, please include the following information: Engle is sold exclusively through authorized Engle Dealers. To find an authorized Engle equipment dealer, call 1-800-547-1906 or email info@engledental.com

- Model and/or Part Number
- Description/Color
- Quantity
- Date Required
- Purchase Order Number
- Shipping Instructions

All prices quoted herein are in US dollars and are subject to change without notice.
All prices listed in this catalog are suggested retail prices.

CUSTOMS

When ordering, you are considered the importer of record and must comply with all laws and regulations of the country in which you are receiving the goods. Contact your local customs office for further information, as customs policies vary from country to country.

RETURNED GOODS

- Before returning the product, call Customer Service for a Return Authorization Number.
Please have the following information when you call:
 - Customer number
 - Part numbers and quantities for all return items
 - The corresponding invoice or order number
- Carefully package all return items to prevent shipping damage.
- Write the Return Authorization Number on the outside of the package and reference the number on any accompanying items.
- Customized products and upholstery returns will not be accepted.
- 20% Restocking Fee


TERMS

Net 30 days on approved accounts; Visa, Mastercard or American Express can be used until credit is approved.

WARNING

Only authorized dealers should attempt to service.
Use of other than authorized dealers will void the warranty

WHAT IS TRAVERSE?



Engle is the only company that makes the "fully programmable 10.25" forward/backward hydraulic Traverse Chair™. The traverse feature allows dentists to position themselves behind the patient comfortably without sacrificing space. Since the seat moves forward as the back of the chair reclines, you can stay put and let the patient come to you.

- FLEXIBLE FOR ANY SIZE OPERATORY
- ENHANCED PATIENT COMFORT FOR LONG PROCEDURES
- ERGONOMICALLY SUITED FOR ANY HEIGHT OF DOCTOR/ASSISTANT
- IDEALLY DESIGNED TO COMBINE WITH REAR DELIVERY FOR 12 O'CLOCK DENTISTRY

engle

**DELIVERY CARTS FOR “LEAN”
AND CABINET FREE OPERATORIES**

STANDARD FEATURES

PLATINUM

- E300 Control Head w/3 Auto HP Holders
- Left/Right Delivery
- StarBoard Work Surface
- 3 Pivot Articulating Control Head
- 3 Pivot Articulating Asst Arm
- Master On/Off Switch
- Asepsis HP Tubing
- Wet/Dry Foot Control
- 3 Way Air-Water Syringe
- H.V.E, Saliva Ejector, Solids Collector
- Individual Water Coolant Systems
- 2L Bottled Water System
- Quick Connect Water Aux
- 4 Storage Bins w/Lid
- Junction Box, 5ft Umbilical
- 5 Year Warranty



Part Number

Part Description

Part Number

Part Description

MOBILE DELIVERY SYSTEMS

ADDITIONAL OPTIONS

P050300	SILVER Mobile Cart DR/ASST
PO50310	GOLD Mobile Cart DR/ASST
P050320	PLATINUM Mobile Cart DR/ASST
PO50330	PLATINUM Mobile Cart ASST Only

P051042	Corian Solid Surface 18"x30"
H014356	6 Pin Handpiece Illumination
H014308	5 Pin Handpiece Illumination
P022755	2nd H.V.E
P070370M	Stainless Steel Tray

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warning.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.



P050135

STANDARD FEATURES

- AS-1 Control Head
- 3-Positions Automatic Handpiece
- Large Work Surface
- Individual Water Coolant Controls
- Asepsis Handpiece Tubing
- Wet/Dry Foot Control
- Handpiece Flush/Purge System
- 3-Way Air-Water Syringe
- Master On/Off Switch
- Solids Collector
- H.V.E.
- Saliva Ejector
- 2 Liter Bottled Water System
- Small Junction Box Air & Water Master Valve/Filter/Regulator Assembly
- Quick Connect Water Auxiliary
- 3-Pivot Articulating Assistant's Arm w/ 4-Position Holder
- Smooth 3-Pivot Articulating Control Head Delivery
- Left/Right Delivery
- 4 Storage Bins w/ Lid
- 5' Long Umbilical
- Assistant's 3-Way Syringe



Part Number	Part Description	Part Number	Part Description
DELIVERY CARTS - DOCTOR'S		DELIVERY SYSTEM OPTIONS	
P050135	Delivery Cart-Doctor & Assistant w/ Table Top	H014356	6 Pin Handpiece Illumination
		H014308	5 Pin Handpiece Illumination
		P099019	Add 4th Automatic Handpiece Position
		P021619	Quick Release For Water Bottle
		P022755	2nd H.V.E.
		P049001	Swivel Tray
		P070370M	Metal Tray Upgrade

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warning.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

P050140



STANDARD FEATURES

- 3-Position Automatic Handpiece Positions
- Individual Water Coolant Controls
- Asepsis Handpiece Tubing
- Wet/Dry Foot Control
- Handpiece Flush/Purge System
- 3-Way Air-Water Syringe
- Master On/Off Switch
- Bottled Water System
- Small Junction Box Air & Water Master Valve/Filter/Regulator Assembly
- Quick Connect Water Auxiliary
- 5' Long Umbilical

Part Number Part Description

DELIVERY CARTS - DOCTOR'S

P050140 Delivery Cart - Doctor's Only

Part Number Part Description

DELIVERY SYSTEM OPTIONS

H014356 6 Pin Handpiece Illumination
 H014308 5 Pin Handpiece Illumination
 P099019 Add 4th Automatic Handpiece Position

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warning.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.



P050145



Shown w/ Optional:
• Metal Tray

STANDARD FEATURES

- Asepsis Handpiece Tubing
- Large Work Surface
- Master On/Off Switch
- Solids Collector
- H.V.E.
- Saliva Ejector
- Bottled Water System
- Small Junction Box Air & Water Master Valve/Filter/Regulator Assembly
- Quick Connect Water Auxiliary
- 3-Pivot Articulating Assistant's Arm w/ 4-Position Holder
- Left/Right Delivery
- 4 Storage Bins w/ Lid
- 5' Long Umbilical
- Assistant's 3-Way Syringe

Part Number

Part Description

DELIVERY CARTS

P050145	Delivery Cart - Assistant's w/ Table Top
---------	--

Part Number

Part Description

DELIVERY SYSTEM OPTIONS

P022755	2nd H.V.E. for Assistant
P070370M	Metal Tray
P049001	Swivel Tray

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warning.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

P050150



STANDARD FEATURES

- 3-Position Automatic Handpiece Holders
- Individual Water Coolant Controls
- Asepsis Handpiece Tubing
- Wet/Dry Foot Control
- Handpiece Flush/Purge System
- 3-Way Air-Water Syringe
- Master On/Off Switch
- Solids Collector
- H.V.E.
- Saliva Ejector
- Bottled Water System
- Small Junction Box Air & Water Master Valve/Filter/Regulator Assembly
- Quick Connect Water Auxiliary
- 5' Long Umbilical

Part Number Part Description

Part Number Part Description

DELIVERY CARTS

DELIVERY SYSTEM OPTIONS

P050150 Delivery Cart - Hygiene

H014356 6 Pin Handpiece Illumination
 H014308 5 Pin Handpiece Illumination
 P099019 4th Automatic Handpiece Position

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warning.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.



STANDARD FEATURES

P071964

- Under Counter Mounted AS-1 Delivery Head
- Left / Right Delivery
- 3 Automatic Handpiece controls
- Individual Water Coolant Controls
- Individual drive air controls
- Asepsis Handpiece Tubing
- Wet/Dry Foot Control
- Handpiece Flush/Purge System
- 3-Way Air/Water Syringe
- Master On/Off Switch
- 4' Umbilical



Part Number

Part Description

DELIVERY SYSTEM

P071964	Rear Delivery, Under Counter Mount - Doctor
---------	---

DELIVERY SYSTEM OPTIONS

H014356	6 Pin Handpiece Illumination
H014308	5 Pin Handpiece Illumination
PO21604	Clean Water Auxiliary W/ Regulator

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

DOCTOR'S & ASSISTANT'S WITH WORK SURFACE



(cabinet is not included in pricing)

PO71438



PO71450

Part Number

Equipment Included

CHOOSE SYSTEM

PO71438	Cabinet Top Mount 12 O'Clock Rear Delivery
PO71450	Cabinet Side Mount 12 O'Clock Rear Delivery

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warning.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.



**A REAL
UN CROWDED
PLEASER**

ENGLE REAR DELIVERY + ENGLE TRAVERSE CHAIR™

Enjoy **UNPARALLELED ERGONOMICS** and **10.25 EXTRA INCHES**
of space behind the patient with the
Engle Traverse Chair and the New Engle Rear Delivery.

For more information, contact 800.547.1906



Part Number **Part Description**

Part Number **Part Description**

DELIVERY SYSTEM

ADDITIONAL OPTIONS

P070936	Engle Dr's Side Delivery System, Wall Mount, E300 Control Head
P070936S	Engle Dr's Side Delivery System, Wall Mount (Short Flex Arm)
P070900	Engle Dr's & Asst Delivery System, Wall Mount AS-1 Control Head
P070902	Engle Dr's Wall Mount Delivery System, AS-1 Control Head
P029355	Engle Panel Mount Control System
P070904	Engle Assistant's Delivery System, Wall Mount

H014356	6 PIN Handpiece Illumination
H014308	5 PIN Handpiece Illumination
P070370M	Stainless Steel Tray
P071655	Wall Mounting Board
P021604	Clean Water Auxiliary System w/Regulator

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warning.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.



STANDARD FEATURES**P070902**

- AS-1 Control Head w/ 3 Auto Handpiece Positions
- Individual Water Coolant Controls
- Asepsis Handpiece Tubing
- Wet/Dry Foot Control
- Handpiece Flush/Purge System
- 3-Way Air-Water Syringe
- Master On/Off Switch
- Small Junction box w/ Air & Water Master Valve /Filter/Regulator
- Smooth 3-Pivot Articulating Control Head Delivery
- Individual Drive Air Adjustments Per Handpiece
- Left/Right Delivery
- 4' Long Umbilical

P070904

- Asepsis Handpiece Tubing
- Solids Collector
- H.V.E.
- Saliva Ejector
- 4-Pivot Articulating Assistant's Arm w/ 4-Position Holder
- Left/Right Delivery
- 4' Long Umbilical
- Assistants 3-way Air/Water Syringe

STANDARD FEATURES**P070936**

- E300 Control Head w/ 3 Auto Handpiece Holders
- Individual Water Coolant Controls
- Master On/Off Switch
- Asepsis Handpiece Tubing
- Handpiece Flush/Purge System
- 3"-Way Air-Water Syringe

P070936S

- E300 Control Head w/ 3 Auto Handpiece Positions
- Short Flex Arm
- Flex Arm w/ Air Brake
- Adjustable Swivel Tray
- Small Junction Box w/ Air & Water Master Valve/Filter/Regulator Assembly
- 8' Long Umbilical
- Wet/Dry Foot Control

STANDARD FEATURES**P070900**

- AS-1 Control Head w/ 3 Auto Handpiece Positions
- Individual Water Coolant Controls
- Asepsis Handpiece Tubing
- Wet/Dry Foot Control
- Handpiece Flush/Purge System
- 3-Way Air-Water Syringe
- Master On/Off Switch
- Solids Collector
- H.V.E.
- Saliva Ejector
- Small Junction box w/ Air & Water Master Valve /Filter/Regulator Assembly
- 4-Pivot Articulating Assistant's
- Arm w/ 4-Position Holder
- Smooth 3-Pivot Articulating Control Head Delivery
- Individual Drive Air Adjustments Per Handpiece
- Left/Right Delivery
- 4' Long Umbilical
- Assistants' 3-way Air/Water Syringe

P029355

- Panel Mount Control w/ 3 Auto Handpiece Holders
- Individual Water Coolant Controls
- Asepsis Handpiece Tubing
- Wet/Dry Foot Control
- Handpiece Flush/Purge System
- 3-Way Air-Water Syringe
- Master On/Off Switch
- Air & Water Master Valve/Filter/Regulator Assembly
- Individual Drive Air Adjustments Per Handpiece
- 6' Long Umbilical



P070937



P070886



P070885

Part Number Part Description

DELIVERY SYSTEM

P070928	Engle 360 Delivery: Radius, Left/Right AS-1 Control Head Delivery System
P070885	Engle 1200 Delivery: Over the Patient, AS-1 Control Head Delivery System w/Cuspidor
P070886	Engle 2200 Delivery: Side Mount, AS-1 Control Head Delivery System
P070912	Engle 2200 Delivery: Chair Mount, AS-1 Control Head, without Assistant's System
P070937	Engle 360 Hygiene Delivery: Radius, Left/Right, AS-1 Control Head w/Vacuum Delivery System (Mounts to Standard Chair Only)

** P070937 NOT COMPATIBLE WITH TRAVERSE CHAIR

STANDARD FEATURES

DELIVERY

- AS-1 Control Head With 2 Automatic Handpieces Saliva Ejector And H.V.E.
- Left/Right Configuration
- Individual Water Coolant
- Master On/Off Switch
- Aspesis Handpiece Tubing
- Handpiece Flush/Purge
- Wet/Dry Foot Control
- 3-way Air/Water Syringe
- Saliva Ejector
- H.V.E.- Solids Collector
- Flex Arm W/ Air Brakes
- 2 Liter Water Bottle With Quick Disconnect
- City Or Water Bottle Selector
- Contoured Box

⚠ WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warning.ca.gov.

⚠ WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.





P069201



P069240 -
shown with optional
P020336 - Asst Syringe
P022755 - 2nd H.V.E



22-507-02

Part Number

Part Description

ASSISTANT'S AND OPTIONS

P069201	Utility Center w/ Assistant - Chair Bracket Needed
P022755	2nd H.V.E. For Assistant
P024343	Engle Chair Bracket
P069240	Utility Center w/ E300 Cuspidor and Assistant- Chair Bracket Needed
P020336	Assistant 3-WAY SYRINGE
P099016	AIR Quick Disconnect
P070904	Assistant's 12 O'clock Delivery System - wall / Cabinet Mounted
P022275	Engle 360 Assistant's Vacuum Delivery System - Left/Right
P069265	Engle 360 E300 Cuspidor Delivery System - Left/Right (Replaces P022275)
22-502-02	Cabinet Vacuum Group - Top Mount
22-507-01	Cabinet Vacuum - Short Arm
22-507-02	Cabinet Vacuum - Long Arm

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warning.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

AVAILABLE AS TRAVERSE OR STANDARD PACKAGE



ASSISTANT'S

PO22275

DELIVERY UNIT

PO70931



#P096265

- Shown w/ Optional:
- Assistant's Touch Pad
 - 2nd H.V.E
 - Bien Air Handpiece
 - Chair Foot Control
 - Metal Tray
 - Twin Stitch Upholstery
 - 360 Cuspidor pictured above

Part Number

Equipment Included

CHOOSE CHAIR

PO96540	Standard Chair
PO96540T	Traverse Chair
PO22275	360 Assistant Arm
PO70931	360 Delivery Unit

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

CHAIR FEATURES – 5 YEAR WARRANTY

- Std. or Traverse Model
- Traverse features 10.25" Movement
- Synchronized back and seat recline
- Dual articulating headrest
- 6 Presets on Touchpad with Entry & Exit
- 60 Degree swivel
- Wideback or Narrowback w/Slings
- Swing down armrest
- Std. Seamless ColorGuard Upholstery
- 120V Chair (220V Upon Request)

LEFT/RIGHT CONVERTIBLE DELIVERY SYSTEM FEATURES – 5 YEAR WARRANTY

- E300 Control Head w/ 3 Auto HP Holders
- Junction Box w/ Air & Water Regulator
- Flex arm with Air Brake
- Master On/Off Switch
- Asepsis HP Tubing
- 2L Water Bottle System w/QD
- 3-Way Air-Water Syringe
- City/Water Bottle Switch
- Handpiece Flush/Purge System
- Quick Connect Water Aux
- Adjustable Swivel Tray
- Individual Water Coolant Controls
- Fully integrate Electric Handpieces, Intraoral Cameras, Scalers, Digital Sensors
- Wet/Dry Foot Control
- 2 FT Umbilical

ASSISTANT FEATURES – 5 YEAR WARRANTY

- Left/Right Convertible
- 3-Way Air-Water Syringe
- H.V.E., Saliva Ejector, Solids Collector
- Upgrade w/ Porcelain Cuspidor #P069265

LED LIGHT – 5 YEAR WARRANTY

- 3-Axis Rotation
- Adjustable LUX (3,000 – 45,000)
- Available Curing Filter
- 5000K Color Temperature

Part Number

Part Description

Part Number

Part Description

TRAVERSE PKG OPTIONS

STANDARD PKG OPTIONS

P096731	Touchpad - Left Side Chair
P096730	Touchpad - Right Side Chair
96-735-01	Touchpad (Control Head Mounted-E300)
96-735-07	Touchpad (Cabinet Mount)
96-735-08	Touchpad (360 Asst, P022275)
P096780T	Chair Foot Control
P099031	Tech Connect USB Upgrade (2 USB Ports & Sensor Holder)
P099033	Audio Jack Integrated into Chair
P096456	I.V. Arm Pad
P099018	Add 4th Handpiece Position
H014356	6 PIN Handpiece Illumination
H014308	5 PIN Handpiece Illumination
P070370M	Stainless Steel Metal Tray
P022755	2nd H.V.E

P096662	Touchpad - Left Side Chair
P096660	Touchpad - Right Side Chair
96-530-01	Touchpad (Control Head Mounted-E300)
96-530-07	Touchpad (Cabinet Mount)
96-530-08	Touchpad (360 Asst, P022275)
P096780	Chair Foot Control (Standard)
P099031	Tech Connect USB Upgrade (2 USB Ports & Sensor Holder) I.V. Arm Pad
P096456	Add 4th Handpiece Position
P099018	6 PIN Handpiece Illumination
H014356	5 PIN Handpiece Illumination
H014308	Stainless Steel Metal Tray
P070370M	2nd H.V.E
P022755	

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warning.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.



SIDE DELIVERY OPERATORY PACKAGE



**DELIVERY UNIT
W/ ASSISTANTS**

PO70925

Part Number

Equipment Included

CHOOSE CHAIR

PO96540	Standard Chair
PO95450T	Traverse Chair
PO70925	Delivery Unit

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warning.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

CHAIR FEATURES – 5 YEAR WARRANTY

- Std. or Traverse Model
- Traverse features 10.25" Movement
- Synchronized back and seat recline
- Dual articulating headrest
- 6 Presets on Touchpad with Entry & Exit
- 60 Degree swivel
- Wideback or Narrowback w/Slings
- Swing down armrest
- Std. Seamless ColorGuard Upholstery
- 120V Chair (220V Upon Request)

DELIVERY SYSTEM FEATURES – 5 YEAR WARRANTY

- E300 Control Head w/ 3 Auto HP Holders
- Utility Center w/quick connect Water Aux
- Flex arm with Air Brake
- Master On/Off Switch
- Asepsis HP Tubing
- 2L Water Bottle System w/QD
- 3-Way Air-Water Syringe
- City/Water Bottle Switch
- Handpiece Flush/Purge System
- Wet/Dry Foot Control
- Adjustable Swivel Tray
- Individual Water Coolant Controls
- 5 FT Umbilical
- Junction Box w/ Air & Water Regulator
- H.V.E, Saliva Ejector, Solids Collector
- Telescopic Assistants Arm w/4 Pos. Holder
- Fully integrate Electric Handpieces, Intraoral Cameras, Scalers, Digital Sensors
- Assistant's 3-way Air-Water Syringe

LED LIGHT – 5 YEAR WARRANTY

- 3-Axis Rotation
- Adjustable LUX (3,000 – 45,000)
- Available Curing Filter
- 5200K Color Temperature

Part Number	Part Description	Part Number	Part Description
TRAVERSE PKG OPTIONS		STANDARD PKG OPTIONS	
P096731	Additional Touchpad - Left Side Chair (Traverse)	P096662	Additional Touchpad - Left Side Chair (Non-Traverse)
P096730	Additional Touchpad - Right Side Chair (Traverse)	P096660	Additional Touchpad - Right Side Chair (Non-Traverse)
96-735-01	Additional Touchpad (Control Head Mounted-E300)	96-530-01	Additional Touchpad (Control Head Mounted-E300)
96-735-07	Additional Touchpad (Cabinet Mount)	96-530-07	Additional Touchpad (Cabinet Mount)
96-735-08	Additional Touchpad (360 Asst, P02275)	96-530-08	Additional Touchpad (360 Asst, P02275)
P096780T	Chair Foot Control (Traverse)	P096780T	Chair Foot Control (Standard)
P099031	Tech Connect USB Upgrade (2 USB Ports & Sensor Holder)	P099031	Tech Connect USB Upgrade (2 USB Ports & Sensor Holder)
P099033	Audio Jack Integrated into Chair	P096456	I.V. Arm Pad
P096456	I.V. Arm Pad	P099018	Add 4th Handpiece Position
P099018	Add 4th Handpiece Position	HO14356	6 PIN Handpiece Illumination
HO14356	6 PIN Handpiece Illumination	HO14308	5 PIN Handpiece Illumination
HO14308	5 PIN Handpiece Illumination	PO70370M	Stainless Steel Metal Tray
PO70370M	Stainless Steel Metal Tray	PO22755	2nd H.V.E
PO22755	2nd H.V.E		

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warning.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.



OVER THE PATIENT WITH CUSPIDOR OPERATORY PACKAGE



DELIVERY UNIT
PO70920

Shown w/ Optional:
• 2nd H.V.E



Shown w/ Optional:
• Metal Tray P070370M



Shown w/ Optional:
• H.V.E

Part Number

Equipment Included

CHOOSE CHAIR

PO96540	Standard Chair
PO95450T	Traverse Chair
PO70920	Cupsidor Delivery Unit

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warning.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.



CHAIR FEATURES – 5 YEAR WARRANTY

- Std. or Traverse Model
- Traverse features 10.25" Movement
- Synchronized back and seat recline
- Dual articulating headrest
- 6 Presets on Touchpad with Entry & Exit
- 60 Degree swivel
- Wideback or Narrowback w/Slings
- Swing down armrest
- Std. Seamless ColorGuard Upholstery
- 120V Chair (220V Upon Request)

DELIVERY SYSTEM FEATURES – 5 YEAR WARRANTY

- E300 Control Head w/ 3 Auto HP Holders
- Utility Center w/quick connect Water Aux
- Porcelain Cuspidor w/ Cup Fill & Bowl Rinse
- Master On/Off Switch
- Flex arm with Air Brake
- 2L Water Bottle System w/QD
- Asepsis HP Tubing
- City/Water Bottle Switch
- 3-Way Air-Water Syringe
- Wet/Dry Foot Control
- Handpiece Flush/Purge System
- Individual Water Coolant Controls
- Adjustable Swivel Tray
- Junction Box w/ Air & Water Regulator
- 5 FT Umbilical
- Telescopic Assistants Arm w/4 Pos. Holder
- H.V.E, Saliva Ejector, Solids Collector
- Assistant's 3-way Air-Water Syringe
- Fully integrate Electric Handpieces, Intraoral Cameras, Scalers, Digital Sensors

LED LIGHT – 5 YEAR WARRANTY

- 3-Axis Rotation
- Adjustable LUX (3,000 – 45,000)
- Available Curing Filter

Part Number	Part Description	Part Number	Part Description
TRAVERSE PKG OPTIONS		STANDARD PKG OPTIONS	
P096731	Additional Touchpad - Left Side Chair (Traverse)	P096662	Additional Touchpad - Left Side Chair (Non-Traverse)
P096730	Additional Touchpad - Right Side Chair (Traverse)	P096660	Additional Touchpad - Right Side Chair (Non-Traverse)
96-735-01	Additional Touchpad (Control Head Mounted-E300)	96-530-01	Additional Touchpad (Control Head Mounted-E300)
96-735-07	Additional Touchpad (Cabinet Mount)	96-530-07	Additional Touchpad (Cabinet Mount)
96-735-08	Additional Touchpad (360 Asst, P022275)	96-530-08	Additional Touchpad (360 Asst, P022275)
P096780T	Chair Foot Control (Traverse)	P096780T	Chair Foot Control (Standard)
P099031	Tech Connect USB Upgrade (2 USB Ports & Sensor Holder)	P099031	Tech Connect USB Upgrade (2 USB Ports & Sensor Holder)
P099033	Audio Jack Integrated into Chair	P096456	I.V. Arm Pad
P096456	I.V. Arm Pad	P099018	Add 4th Handpiece Position
P099018	Add 4th Handpiece Position	H014356	6 PIN Handpiece Illumination
H014356	6 PIN Handpiece Illumination	H014308	5 PIN Handpiece Illumination
H014308	5 PIN Handpiece Illumination	P070370M	Stainless Steel Metal Tray
P070370M	Stainless Steel Metal Tray	P022755	2nd H.V.E
P022755	2nd H.V.E		

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warning.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.





STANDARD FEATURES

- 3 Axis Rotation Light Head
- Focal Point: 3"x 6" @ 27"
- Adjustable LUX (3,000 - 45,000)
- No Infrared Light
- 5000K Color Temperature
- 10+ Years of LED Life
- Low Voltage
- Available Curing Filter

Part Number

Part Description

LED LIGHTS

PO24730	360 MAIA LED LIGHT (RADIUS), Must Be Ordered With Engle Chair
PO24720	320/310 MAIA LED Light, Utility Center Mounted
PO24710	MAIA LED Post Mounted Light (Need Chair Bracket & Post)
PO24145	32" Engle Light Post (P024710)
PO24140	40" Engle Light Post (P024710)
PO24681	Antipolimerization Filter (Curing Filter)

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warning.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.


STANDARD FEATURES

- HDMI Cable Factory Installed
- Mounts to Engle Chair or Utility Center
- VESA Mount Compatible, size 100mm and 75mm
- Articulating Monitor Mount with Quick Click Positioning
- Monitor Recommendations:
 - 27" Max Monitor Size
 - DC Power Only

** P024620:
Monitor Not Included



P099032



P024610

Part Number

Part Description

MONITOR MOUNTS

P024610	360 Monitor Mount (Radius)
P099032	320/310 Utility Center Mounted Monitor

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warning.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

ENGLE 300 & ENGLE 300 TRAVERSE



525LB LIFT CAPACITY -
ALL HYDRAULIC MOVEMENT



P096540 - Standard Chair
P096540T - Traverse Chair

10.25" TRAVERSE FORWARD/BACK

Part Number

Part Description

CHOOSE CHAIR

P096540T	Engle Traverse Hydraulic Chair
P096540	Engle Standard Hydraulic Chair

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

STANDARD FEATURES

STANDARD CHAIR

- Thin Back w/ Lumbar Support
- Wide Back OR Narrowback w/ Slings
- Seamless ColorGuard Upholstery
- Anatomically Correct Back Pivot Point
- 6 Programmable Positions/4 Manual
- Touchpad (Specify Location - Lt or Rt Side)
- Synchronized Back & Toe Board Recline
- Smooth, Quiet All Hydraulic System
- 60° Seat Rotation
- Quick Release Double Articulating Headrest
- Left/Right Swing Down Upholstered Arms
- Programmable Exit

TRAVERSE CHAIR

- Thin Back w/ Lumbar Support
- Wide Back OR Narrowback w/ Slings
- Seamless ColorGuard Upholstery
- Anatomically Correct Back Pivot Point
- 6 Programmable Positions/6 Manual
- Touchpad (Specify Location - Lt or Rt Side)
- Synchronized Back & Toe Board Recline
- Smooth, Quiet All Hydraulic System
- 60° Seat Rotation
- Quick Release Double Articulating Headrest
- Left/Right Swing Down Upholstered Arms
- Programmable Exit
- 10.25" Hydraulic Traverse Movement

Part Number	Part Description	Part Number	Part Description
TRAVERSE CHAIR OPTIONS		STANDARD CHAIR OPTIONS	
P096731	Additional Touchpad - Left Side Chair Mount (Traverse)	P096662	Additional Touchpad - Left Side Chair Mount (Non-Traverse)
P096730	Additional Touchpad - Right Side Chair Mount (Traverse)	P096660	Additional Touchpad - Right Side Chair Mount (Non-Traverse)
96-735-01	Additional Touchpad (Control Head Mounted-E300)	96-530-01	Additional Touchpad (Control Head Mounted-E300)
96-735-07	Additional Touchpad (Cabinet Mount)	96-530-07	Additional Touchpad (Cabinet Mount)
P096780T	Chair Foot Control (Traverse)	P096780	Chair Foot Control (Standard)
P099031	Tech Connect USB Upgrade (2 USB Ports & Sensor Holder)	P099031	Tech Connect USB Upgrade (2 USB Ports & Sensor Holder)
P099033	Audio Jack Integrated into Chair	P099033	Audio Jack Integrated into Chair
P096456	I.V. Arm Pad	P096456	I.V. Arm Pad
P099018	Add 4th Handpiece Position	P099018	Add 4th Handpiece Position
H014356	6 PIN Handpiece Illumination	H014356	6 PIN Handpiece Illumination
H014309	5 PIN Handpiece Illumination	H014309	5 PIN Handpiece Illumination
PO70370M	Stainless Steel Metal Tray	PO70370M	Stainless Steel Metal Tray
PO22755	2nd H.V.E	PO22755	2nd H.V.E

 **WARNING:** This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

 **WARNING:** This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

STANDARD ULTRALEATHER COLOR OPTIONS

UltraLeather is the ultimate in softness and durability. Performance fabric redefined.


Engle offers 35 Standard Colors along with the full line (100+) of "Custom ULTRALEATHER" colors.



STANDARD COLORGUARD OPTIONS

ColorGuard features a durable, 36 oz. vinyl construction that is resistant to microbial and bacterial attack. Featuring the industry-leading PreFixx protective finish.

Engle offers 10 Standard Colors along with the full line (70+) of "Custom ColorGuard" colors.



Request Free Upholstery Samples Online at www.engledental.com or email info@engledental.com

Please note printed catalog and computer images of color shades are approximate only, actual material swatch will vary in shade/tone from colors shown.



Twin Switch



Summit



Seamless

Part Number

Part Description

CHOOSE UPHOLSTERY: OPTION OF WIDEBACK OR NARROWBACK W/ SLINGS*

*NARROWBACK AVAILABLE WITHOUT SLINGS UPON REQUEST

P096683	Wideback - Seamless ColorGuard Upholstery
P096685	Wideback - Upgrade to Seamless Ultraleather Upholstery
P096687	Wideback - Upgrade to Twin Stitch Ultraleather Upholstery
P096824	Wideback - Upgrade to Summit Ultraleather Upholstery
P096685HM	Wideback - Upgrade to Seamless Heat & Massage Ultraleather Upholstery
P096687HM	Wideback - Upgrade to Twin Stitch Heat & Massage Ultraleather Upholstery
P096825HM	Wideback - Upgrade to Summit Heat & Massage Ultraleather Upholstery
P096684	Narrowback w/Slings - Seamless ColorGuard Upholstery
P096686	Narrowback w/Slings - Upgrade to Seamless Ultraleather Upholstery
P096688	Narrowback w/Slings - Upgrade to Twin Stitch Ultraleather Upholstery
P096826	Narrowback w/Slings - Upgrade to Summit Ultraleather Upholstery
P096686HM	Narrowback w/Slings - Upgrade to Seamless Heat & Massage Ultraleather Upholstery
P096688HM	Narrowback w/Slings - Upgrade to Twin Stitch Heat & Massage Ultraleather Upholstery
P096827HM	Narrowback w/Slings - Upgrade to Summit Heat & Massage Ultraleather Upholstery
EMBR	Custom Embroidery - Image file Required with Order
DUAL	Dual Tone or Custom Upholstery

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warning.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.



STANDARD - P097092

- Comfortable 4" Thick, 15" Diameter Seat Cushion
- Many Options Allow Ability to Customize
- Contour Backrest (Optional)
- Various Height Options
- Swing Arm or Ratched Body Support



CONTOUR - P097094

- Contoured 18" x 18" Seat Cushion
- 24" Powder-Coated Cast Aluminum Base
- Dual Adjustable Lumbar Backrest Adds Comfort
- Modern Design Fits Any Office



ERGONOMIC - P097096

- Accommodates Front to Back Reaching without Sacrificing Lumbar Support
- Exclusive Use of HybriGel Foam
- Independently Adjust Backrest to All Greater Range of Movement
- Provides Correct Weight Distribution



STANDARD - P097093

- Comfortable 4" Thick, 15" Diameter Seat Cushion
- Many Options Allow Ability to Customize
- Contour Backrest (Optional)
- Various Height Options
- Swing Arm or Ratched Body Support



CONTOUR - P097095

- Contoured 18" x 18" Seat Cushion
- 24" Powder-Coated Cast Aluminum Base
- Dual Adjustable Lumbar Backrest Adds Comfort
- Modern Design Fits Any Office



ERGONOMIC - P097097

- Accommodates Front to Back Reaching without Sacrificing Lumbar Support
- Exclusive Use of HybriGel Foam
- Independently Adjust Backrest to All Greater Range of Movement
- Provides Correct Weight Distribution

Part Number

Part Description

DOCTOR & ASSISTANT STOOLS

P097092	Doctor Stool - Standard, Vinyl Upholstery
P097093	Assistant Stool- Standard, Vinyl Upholstery
P097094	Contour Doctor Stool - Vinyl Upholstery
P097095	Contour Assistant Stool- Vinyl Upholstery
P097096	Ergonomic Doctor Stool- Vinyl Upholstery
P097097	Ergonomic Assistant Stool - Vinyl Upholstery
P097092UL	Doctor Stool - Standard, Ultraleather Upholstery
P097093UL	Assistant Stool- Standard, Ultraleather Upholstery
P097094UL	Contour Doctor Stool - Ultraleather Upholstery
P097095UL	Contour Assistant Stool- Ultraleather Upholstery
P097096UL	Ergonomic Doctor Stool- Ultraleather Upholstery
P097097UL	Ergonomic Assistant Stool - Ultraleather Upholstery

WARNING: This product can expose you to chemicals including Titanium Dioxide, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warning.ca.gov.

WARNING: This product can expose you to Lead and Lead Compounds, which is known to the State of California to cause cancer and/or birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.






ENGINEERED TO PERFORM



SERVICE & SUPPORT





Handcrafted in Hillsboro, Oregon, Engle products are engineered and manufactured under our own roof.





CUSTOMIZABLE

With access to a team of technicians and the manufacturing team itself, Engle provides fast and efficient problem-solving.


With your choice of custom colors, fabrics, and upholstery patterns, you can design your chair to suit your style.





P070931





P070350





P024610




P069265




P022275





P070920





310 UNIT ARM IDENTICAL TO 320




P070925




P070935




P070900 - DR'S D.P. WALL MOUNT




P070902 - DR'S WALL MOUNT



P070904 - ASSISTANT'S WALL MOUNT



P071964 - CABINET MOUNTED UNIT



AS SEEN IN:

DENTAL PRODUCT SHOPPER

ENGLE DENTAL SYSTEMS

Delivering a Cutting-Edge Operator Solution

Engle's range of dental chairs, lights, delivery systems, and supporting equipment were designed to be used together to create the complete, ergonomic dental operator

Whether you're remodeling an operator or building a new office, choosing the right equipment can be tricky. There are dozens of factors to consider such as how much space is available; where is the best location for instruments, lighting, and monitors; and whether you practice left-, right-, or four-handed dentistry. Plus, there's the consideration of how the setup affects patient comfort, treatment

time, workflow, and the overall team experience.

No matter what space and conditions your practice has to work with, Engle Dental Systems can provide customized, quality operator equipment that meets your unique needs. For the past 50 years, Engle has developed dental chairs, delivery units, and other equipment with a focus on state-of-the-art quality and ergonomics. All equipment is designed and manufactured at its Hillsboro, OR, headquarters, and quality control is handled in-house.



CONNECT TO THE FUTURE

Engle always has the future of digital dentistry in mind. The **TechConnect Operator Package** is a first-of-its-kind system that powers a practice's existing digital technology through prewired HDMI cables, 2 USB ports for internal and external use, and a digital sensor holder. This offers clinicians simple chairside access to tools and technology like digital sensors, handpieces, intraoral cameras, and monitors.

Optimize the Operator

Remodeling or designing a new room begins with the dental chair. Engle's fully programmable, 10" forward/backward Traverse Chair allows clinicians to position themselves behind the patient comfortably without sacrificing space. Built with advanced hydraulic controls, the seat moves forward as the back of the chair reclines so dentists can let the chair position itself. "The chair's ability to traverse was huge for our practice," said hygienist and DPS evaluator Debbie Alo, who found the Traverse Chair to be an ideal solution for her small hygiene room.

In addition to the dentist's choice of color, fabric, and upholstery pattern, the Standard and Traverse chair models feature lumbar support, a double articulating headrest, and other features to help patients relax. To complement the chair, Engle carries four operator LED light types and a wide variety of delivery units and systems, with options for side delivery, wall-mounted, and carts. Engle customers also receive access to a world-class team of technicians and manufacturers for fast and efficient problem-solving.



www.engledental.com

Handcrafted in Oregon



USA

Engle Dental Systems LLC
7205 NE Evergreen Parkway, Suite 100
Hillsboro, Oregon 97124
Tel: (800) 547-1906 • Fax: (503) 357-7286
info@engledental.com

P242000 6.0