


# QuickRay

Direct USB



## Features & Highlights:

- CMOS technology with scintillator
- 2 Sensor Sizes
- 4096 gray scale levels (12 bit image)
- Direct USB connection
- No data box
- No external power connection
- Built in USB 2.0 Interface
- USB 2.0 Plug and Play
- Two sensors can be operated
- Image transfer and control via the PC
- 3 meter USB sensor cable
- 12 bit image digitalization
- 16 bit uncompressed image file save
- Well defined software interface Standard or 2x2 binned imaging modes
- Intelligent Auto Trigger

## Technical Specifications:

- Theoretical resolution 26.3 lp/mm
- Real resolution >20 lp/mm or 0.05 mm image detail
- Signal to noise ratio: 40 dB
- Scintillator: CSI depositing on fiber-less sub-strate

**e2v**

“Inside” Components

## Most Software Integration on the Market!

Dexis, Gendex-VixWin, Carestream-Kodak, Schick-Eaglesoft, Apteryx & MacPractice!



**Special Price!**

**Size #2 Sensor \$3995**

**Size #1 Sensor \$3550**

Please Contact Randy Ahlquist @ Integrity Dental (303) 888-5500