

vakker[®]

2025-2026



vakker[®] Dental Inc



: www.vakkerdental.com



: info@vakkerdental.com



: (888) 516 0083



1. Digital Radiography	<1-3>
2. Endodontics	<4-12>
3. Ultrasonic Scaler & Tips	<13-15>
4. Air Polisher	<16-18>
5. Curing Light	<19-20>
6. Handpiece	<21-23>
7. Sterilization Pouches	<24-25>
8. Dental Electric Motor & Handpiece	<26-28>
9. Diode Laser	<29-30>
10. Woodpecker Implant Motor	<31-32>
11. Ultrasonic Bone Surgery	<33-34>
12. Portable Micromotor	<35-36>
13. Automatic Crown Remover	<37-38>
14. Dental Anesthesia Delivery Device	<39-40>
15. Dental Microscope	<41-42>
16. Teeth Whitening	<43-55>

01

Digital Radiography

Woodpecker Dental Ai Ray Pro & Ai Ray Lite Machine

■ Ai Ray Pro 27-98

Japanese Canon made 0.4mm X-Ray tube, high frequency, more stability, safe photographing, more detailed clearer image and more reliable.



■ Ai Ray Lite 27-116

Japanese Canon made 0.4mm X-Ray tube, high frequency, more stability, safe photographing, more detailed clearer image and more reliable.



vakker® **Dental Practice Helper!**

i-Sensor

Clearer & More durable



H1 (38.5x25x4.5mm/1.52x0.98x0.18inch)

H1.5 (43.5x29.2x4.5mm/1.71x1.15x0.18inch)

H2 (45x31.6x4.5mm/1.77x1.24x0.18inch)

i-Scan 27-47

Woodpecker i-Scan
Imaging Plate Scanner



2 year warranty Free Software Twain Compatible Free Technical Support

Twain Driver Compatible with Apteryx DTX Studio, Apteryx XrayVision, Apteryx XVCapture, CliniViewXV, ClioSoft/SotalImage, CurveDental/Curve Hero, Dentrax Ascend, DentiMax, Dexis, EVAsoft, Harmony Jazz, iDentalSoft, i-Dixel, MiPACS, OpenDental, Oryx Imaging, Planmeca Romexis, SoPro Imaging, Sota Cloud, Sota Imaging, Tab32, Teledentix, TigerView, Triana, Visix Imaging, XDR etc

02

Endodontics


vakker® Dental Practice Helper!

Woodpecker Ai-Motor Endo Motor with Apex Locator



Ai-Motor 27-60

Design by Dr. Yoshi Terauchi



T-Mode Safe Ledge Bypassing



Woodpecker Worldwide Unique Endo Motor with Vibration

812 MT 27-130

Endo Motor



vakker® Dental Practice Helper!

Woodpecker Endomatic & Endo Free Endo Motor with Apex Locator

EndoMatic 27-2-RCM

Specification:

Speed range : 100-1000 RPM

Torque range : 0.4-5.0 N.cm

Integrated length determination
The EndoMatic combines with the length measuring function.



Endo Free 27-52

- Lightweight handpiece
Easy preparation
- Brushless motor
Strong power
- 10 times the service life of brush motor



VK88 9 Program Reciprocating LED Endo Motor



vakker®

VK88

Specification:

1. Electric motor
2. Drive-16:1 electric contra angle is available
3. Easy operating with wide LCD screen
4. 9 program setting can be memorized, adjustable speed and torque
5. Auto reverse to protect the file
6. Auto power off and memory function
7. 1200mAh high capacity Li-ion chargeable battery



vakker[®] Dental Practice Helper!

Ai-Pex Endo Apex Locator w/ Pulp Tester & DPEX V Apex Locator

Ai-Pex 27-51

APEX LOCATOR WITH PULP TESTER



DPex V 27-48

Apex Locator



Woodpecker Fi-G/Fi-G/Fi-E Gutta-Percha Obturation System



Fi-G 27-49-G

Fi-P 27-49-P

Fi-E 27-49-E

GUTTA-PERCHA OBTURATION SYSTEM



Woodpecker Endo 3 Ultrasonic Endo Activate Device

Endo 3 27-74


- LCD screen + Ring switch button
- Dual-mode setting
- Three power levels
- Four time modes
- Wireless control, handy to operate




Woodpecker Gutta-percha Obturation Pen & Gun




Fi-E Tips



Gun Tips



Pen Tips



Woodpecker R1 Plus Endo Ruler

R1 Plus 27-49-R1-Plus

Gutta-percha Cutter

Designed by international dental experts
 Conducting 220 precise manufacturing processes
 on the equipment imported from Germany and Japan



Vakker Plugger

PLUGGER 27-49-FI-N



- 4 models with different tip diameters
 35-70
 40-80
 50-100
 60-120

- The smaller end is NiTi alloy tip and the other end is stainless steel tip.
 Good rust resistance.
 Super elastic NiTi alloy can easily restore to its original state after 500 times of bending.
 Two laser marking rings of 10 mm and 15 mm engraved on the NiTi tip enable easier identification of the condensation depth.

- 78 precise machining processes are conducted on each plugger with CNC machining equipment imported from Germany and Japan

- Taper accuracy up to 0.02mm



- Ergonomically designed handpiece allows easy operation and comfortable holding



Vakker® Endodontic Files

IV ONE GOLD

(Reciprocating Files)

Description

Can Replace the WaveOne Gold Reciprocating Files

Size: Small 20.07 Primary 25.07 Medium 35.06 Large 45.05

Length: 21mm 25mm 31mm

Package: 3/pk single package; 4/pk assorted package



V-Taper Gold

(NiTi Rotary Files, Protaper Gold Users)

Packing : 6pcs/pack single size / assorted

Size : Shaping Files: #SX, #S1, #S2

Finishing Files : #F1, #F2, #F3, #F4, #F5

Length : 21,25,31mm

Taper : Variable taper

Cross Section : 

Features

New generation of raw material; Advanced technical process

Shocked Toughness

200%-500% Fracture Resistance improved

SX 0.19mm S1 0.18mm S2 0.20mm

F1 0.20mm F2 0.25mm F3 0.30mm F4 0.4mm F5 0.5mm

Taper:04-08

Torque:5.0 N/CM(SX,S1); 1.5 N/CM(S2,F1); 3.0 N/CM(F2,F3,F4,F5)

Recommended Rotary Speed: (260-350 rpm)

Recommended Temp: ≤126



Vakker® Endodontic Files

V BLUE

(Vortex Blue Users)

Q: Which system is advantageous in curved canals ?

A: V Blue !

15	15/.04	15/.06	20	20/.04	20/.06
25	25/.04	25/.06	30	30/.04	30/.06
35	35/.04	35/.06	40	40/.04	40/.06
45	45/.04	45/.06	50	50/.04	50/.06

Material: NiTi **Standard:** Blue **Length:** 21,25,31 mm

Rotary speed: 250-360 rpm **Torque:** 2.0~3.0 N/cm

Package: 6 pcs per box/blister assorted/single size



K FILES

006	008	010	015	020	025	030	035	040	045	050	055	060	070	080
-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----	-----

Material: Stainless Steel/NiTi **Package:** 6 pcs per box assorted/ single size

Size: ISO sizes **Sterilization:** ≤126°C **Length:** 21,25,28,31 mm



Hand use



Engine-driven



C FILES


006	008	010	015	020
-----	-----	-----	-----	-----

Material: Stainless Steel **Package:** 6 pcs per box assorted/ single size

Size: ISO sizes **Sterilization:** ≤126°C **Length:** 18,21,25 mm

Features and Benefits:

- Ensures the patient's safety
- Reduces the risk of instrument breakage
- Twisted file design for greater strength



Hand use




Vakker® Endodontic Files

BARBED BROACHES



Material : Stainless Steel **Package :** 10 pcs per blister assorted/ single size

Size : ISO sizes **Sterilization :** ≤126°C **Length :** 21,25,28 mm



Hand use



PEESO REAMERS



Material : Stainless Steel **Length :** 28,32 mm

Max. Rotary speed : 1200 min⁻¹ **Package :** 6 pcs per box



GATES DRILLS



Material : Stainless Steel **Length :** 28,32 mm

Max. Rotary speed : 800 min⁻¹ **Package :** 6 pcs per box



- Advanced material: differed from 304 stainless steel, imported from South Korea
- Higher hardness martensitic grade for cutting instrumentation
- Higher fracture resistance compared with other brands.



03

Ultrasonic Scaler & Tips

DTE D600 LED Ultrasonic Scaler



D600 LED 27-19-5

Ultrasonic Endodontic Expert
 Intelligent power output enables
 painless treatment--Sine wave drive;
 Intelligent degenerative feedback;
 Smart touching system--Ultra-sensitive
 & Waterproof design;

Features:

- 1 · Sealed handpiece autoclavable.
- 2 · Scaling, Perio, Endodontic
- 3 · Automatic frequency tracking.



DTE D5 LED Ultrasonic Scaler

D5 LED 27-19-2

Feature:

Three functions: Scaling, Perio, Endo.

Technical Specification:

Power Input: 110V~50Hz/60Hz 1.2A (max)

Water Pressure: 0.01MPa~0.5MPa

Main Unit Input: 24V~1.3A

Dimension: 188mm*159mm*81mm

Output Power: 3W~20W

Frequency: 28kHz±3kHz

Main Unit Weight: 0.65kg

Adapter Weight: 1.2kg

1 DTE D5 LED Ultrasonic Scaler 6 tips

(2*GD1T,1*GD2T,1*GD4T,1*PD1T,1*ED1T)



DTE S6 LED Ultrasonic Scaler

S6 LED 27-50-S6

Main unit input: 24v 50hz 110 v
 Water pressure: 0.01MPa-05MPa
 Frequency: 28khz + 3khz
 Main unit weight: 650g
 Dimensions: 197mmx174mmx74mm
 S6 LED Ultrasonic Scaler with 7 Tips
 (GD5, GD12, PD4, PD7, PD12, ED60*2
 Tips)



DTE D7 LED Ultrasonic Scaler

D7 LED 27-19-4

The ultrasonic scaler is composed of the electro circuit, water spray, and ultrasonic transducer.

Application:

Woodpecker DTE D7 LED ultrasonic scaler is used for dental calculus elimination and root canal treatment.

Functions:

Woodpecker D7 LED ultrasonic scaler has scaling, period endo, and auto-water supply functions.

Technical Specifications:

Power input: 220 - 240V~ 50Hz/60Hz / 110V~ 50Hz/60Hz
 Main unit input: 24V~ 50Hz/60Hz 1.3A
 Output power: 3W~20W
 Frequency: 28kHz±3kHz
 Main unit weight: 1.35kg
 Adapter weight: 1.2kg
 Dimensions: 237mm × 206mm × 95mm
 D7 LED Ultrasonic Scaler with 8 Tips (GD1*2,
 GD4, GD5 PD1,PD3, ED1,PD4D Tips)



04

Air Polisher

Woodpecker AP-A Dental Air Polisher & Scaler

AP-A 27-114

Main technical specification:

Power adapter input: 100-240V~ 50/60Hz 1.1A
 Model: USE48- 300130SPA3
 Main unit input: 30V 1.3A
 Main unit insurance: T1.6AL 250V
 Air inlet pressure: 5.5bar ~ 7.5bar(0.55MPa ~ 0.75MPa)
 Main unit weight: 2.0kg
 Main unit size: length×width×height 270mm×170mm×90mm

Faster and more convenient scaling
 Two-in-one: Supragingival & subgingival mode
 Ergonomic design
 Fully transparent funnel-shaped powder tank
 Compact and lightweight



Woodpecker Dental AP-H L & AP-H Plus L Air Polisher

AP-H Plus L

AP-H L

Main technical specification:

Input water pressure: 0.7bar to 2.2bar (70-220kPa)
 Input air pressure: 3.5bar to 4.5bar (350-450kPa)
 Working mode: continuous operation, Intermittent operation
 Weight: About 0.150kg
 Environment temperature: + 10℃ ~ +40℃
 Relative humidity: 30% ~ 75%
 Atmospheric pressure: 50kPa ~ 106kPa



AP-H L: 1pc Supragingival Handpiece
 AP-H Plus L=AP-H L+AP-2 Plus L Subgingival Handpiece+ 10pc Subgingival Nozzles



PT-B

Five - in - One

Comfort Reigns Convenience Rules

- Air polishing
- Periodontal treatment
- Endodontic treatment
- Restoration preparation
- Implant maintenance

27-112 PT-B

State-of-the-Art Perio Treatment Protocol Unique Guided Biofilm Therapy (GBT)

Air polishing

- Realize both supra-and subgingival scaling
- 12 levels of water volume and power adjustment for precise treatment
- 4 temperature settings for added comfort

Ultrasonic scaling

- Three modules: G P E
- Four functions:

Periodontal treatment	Endodontic treatment
Cavity preparation	Implant maintenance
- To offer a one-stop ultrasonic solution

System setting

- Central processing interface
- Multiple mode selection
- One-button switch



05

Curing Light

vakker® Dental Practice Helper!

Woodpecker O-Star & O-Light 1 Sec Curing Light

1 year warranty

■ O-star 27-54

■ O-star Pro 27-54

■ O-Light 27-4-14

■ O-Light Max 27-4-15

O-Star/O-Star Pro Wide Spectrum
Is Curing Light from Woodpecker



VK-018 1 Sec Broad Spectrum Curing Light

■ VK-018

Key Features:

- ◆ Light Intensity: Upto 2300 mW/cm²
- ◆ Wave Length: 370-550 nm

Cure 2mm resin within 1 second

Mode:

- P1 Turbo 1s 3s;
- P2 Normal 5s 10s 15s 20s

Luxury Metal Tips Anti Fall Anti Cracking &
Fast Heat Dissipation
Continue using upto 160 hours

Color available

- White (order#27-4-4-WT)
- Black (order# 27-4-4-BK)

1 year warranty



06

Handpiece

Vakker High Speed Air Turbine Handpiece

■ V1

V1 Plus Non-Optic High Speed Handpiece

Japanese Turbine

6 months warranty



VK200/VK-45/VK-PM01 Handpiece

More Powerful, More Durable, Less Cost

12 months warranty



works with
4 & 6 Hole



■ VK-200
(with fiber optic)

■ VK-45

■ VK-PM01
(mini head)

■ VK100

VK100 High Speed Turbine E-generator Handpiece

Fiber Optic German Turbine

12 months warranty



■ VCQC01

Vakker Quick Coupler

Built-in Generator Works for 4 & 6 Holes Compatible for Kavo



vakker[®] Dental Practice Helper!

■ VKST16

E type External Water Spray Low Speed Handpiece & Air Motor Set

E-type 16:1 Straight Nose Cone Reduction Speed Handpiece



1 Year Warranty



1 Year Warranty

Dental Hygiene 4 Hole Prophy Air Motor & Handpiece Set/Prophy Angles

■ VK-PP02 ■ VK-PP01

■ Blue Firm Prophy Angle ■ Purple Soft Prophy Angle

■ Prophy Mini

12 months warranty



07

Sterilization Pouches

Vakker Sterilization Pouches



Sterilization Pouches

2 1/4" x 5" (2 1/4" x 4") (57x100mm)



Sterilization Pouches

3 1/2" x 10" (3 1/2" x 9") (90x230mm)



Sterilization Pouches

2 3/4" x 10" (2 3/4" x 9") (70x230mm)



Sterilization Pouches

4 1/4" x 12" (4 1/4" x 11") (108x280mm)




Sterilization Pouches

3 1/4" x 12" (3 1/4" x 11") (83x280mm)




Sterilization Pouches

5 1/4" x 11" (5 1/4" x 10") (135x260mm)



Sterilization Pouches

3 1/2" x 6 1/4" (3 1/2" x 5 1/4") (90x135mm)



Sterilization Pouches

7 1/2" x 14" (7 1/2" x 13") (190x330mm)

08

Dental Electric Motor & Handpiece

vakker® Dental Practice Helper!

Woodpecker® ES5 Electric Surgical Motor w/ 1:5, 1:3 handpiece

ES5 27-109

5.0-inch LCD
color touch
screen

2
modes

8
speed
ratios

6
preset
programs

Automatic saving parameters

**2 year on main unit & 1 year
on handpiece Warranty**



Woodpecker Electric Handpiece with Fiber Optic

1 Year Warranty



1:4.2

WJ-45L



1:3

WJ-SG45CL



1:5

WJ-45TL

vakker[®] Dental Practice Helper!

Woodpecker MT2 & MT3 Dental Electric Motor

2 Year Warranty



■ MT3 27-99



■ MT2 27-108



MT3 Working independently, with air compressor system, come with foot pedal, no need to connect with dental unit. so it can be used without an air compressor/ delivery unit connection.

Vakker Dental Electric Motor



■ VKEM01

- Faster Tooth Preparation
- Less Splatter
- More Stable Power Supply
- 2 Year Warranty



■ WJ-45L
1:4.2 Handpiece



■ WJ-45TL
1:5 Handpiece

09

Diode Laser

Woodpecker Dental Smart Diode Laser



Smart Laser 27-66SM

Innovative Two/Three Wavelengths Blue Semiconductor Laser Therapy Device

Can be used for the fields of dental soft tissue surgery, endodontic sterilization, periodontal sterilization, peri-implantitis, low intensity laser therapy, oral ulcer, and teeth whitening.

Remarkable treatment effect: There is basically no bleeding during operation, maximally minimizing the postoperative swelling, reducing pain of patients, and accelerating the recovery.

Power adapter input: 100-240Vac, 50/60Hz, 2.5A

Main unit input: 15V 6.0A

Wavelength and power:

a) 976± 20nm : 0.2W-4W(CW), Peak Power 5W

chopped 1Hz to 20KHz

chopped mode:5us – 0.9S

b) 650± 20nm : 25mW-200mW(CW)

c) 450± 20nm: 0.2W-3W(CW)

chopped 1Hz to 20KHz

chopped mode:5us – 0.9S

Laser classification:

a) 976 nm: Class 4

b) 650 nm: Class 2

c) 450 nm: Class 4

Aiming beam: 650 ± 20 nm / Pmax<5mW (class 1)

Aiming beam: 650 ± 20 nm /Pmax<5mW

Divergence:12.7°(According to IEC 60825-1: 2014)

Rechargeable battery: 11.1V/2600mAh x2 (57.7Wh)

Time consumption for charging: about 4h (5 hours for first charging)

Diverging half angle: 0.22 rad/ 12.6 °



Before treatment



Under treatment



Four weeks later

10

Woodpecker Implant Motor

Woodpecker Implanter Plus

■ Implanter Plus 27-70

Specification:

Torque Range/Gear20:1: 5N.cm-80Ncm

Device for intermittent operation: 3 min ON,10 min OFF

Power supply voltage: 200-240V~

Power supply frequency: 50/60Hz

Software version: Implant-V1

Input power: 150VA

Fuses: 2×T1.6AL 250V

Applied parts: Contra-angle handpiece

Maximum Temperature: 41.8℃

Maximum water volume: 135ml/min

2 Year warranty on main unit & 1 Year warranty on handpiece



Woodpecker Implant Air w/ 20:1 Fiber Optic Handpiece

■ Implant Air 27-139

2 Year warranty on main unit &

1 Year warranty on handpiece

The Woodpecker Implant Air is the latest and most advanced version of implant motor from Woodpecker which is completely redesigned to be simple, precise and reliable for your implant surgery. Implanter Air will be your favorite assistant of high-precision dental implanting treatment.



11

Ultrasonic Bone Surgery

Woodpecker US-II LED Ultrasonic Bone Surgery

■ US-II LED 27-61

Specification:

1 x Main Unit,
 2 x LED Handpieces with cord,
 14 x Tips,
 Basic Kit - US1, US2, US3, US4, US5
 Sinus Lift Kit - UL1, UL2, UL4, UL5
 Osteotomy Kit - US1, US2, US1L, US1R
 Implant Prep Kit - UI2, UI3, UL2

2 x Tip Holders,
 1 x Torque Wrench,
 1 x Foot Switch,
 1 x Pump connector,
 1 x Surgical Tray,
 1 x Suit case,
 1 x Rod,
 1 x connector etc.



Woodpecker Surgic Smart Ultrasonic Bone Surgery

■ Surgic Smart 27-67

Specification:

Device in accordance with Directive 93/42EEC.
 According to EN60529: IPX1 (device) IPX6 (foot pedal)
 Software version: BoneSurgical4-V1.0.0
 Device for intermittent operation: 60s ON, 10s OFF
 Power-supply voltage: ~100V-240V 50Hz/60Hz 170VA
 Fuses: 2x1.6AT 250V
 Working frequency: 24kHz~36kHz
 Flow: 25~110ml/m
 Applied parts: Tips



12

Portable Micromotor

vakker[®] Dental Practice Helper!

Vakker Portable Brushless Micromotor

Fit for dental, manicure, podiatry, jewelry, handicraft processing and so on.



VKMM1



VKMM1E

Speed: 1,000-35,000 RPM Torque: up to 9.6 N.cm

VK-MM1 Power Control Box		VK-MM1 Charger	
Maximum Current	Input voltage: DC 12.6V 1A Output voltage: 12.6V 4A	Maximum Current	Input voltage: 100-240V 47-63HZ Output voltage: 12.6V 1A
Weight	322g		
Dimensions	132*65*61mm		
Battery capacity	2600mah		

VK-MM1 Handpiece	
Maximum speed	35000rpm
Voltage	12.6V
Weight	162g
Dimensions (mm)	L140 D16 H18



VK-MM1E Handpiece	
Maximum speed	35000rpm
Voltage	12.6V
Weight	120g
Dimensions (mm)	L140 D16 H18


13

Automatic Crown Remover

vakker[®] Dental Practice Helper!

Automatic Crown Remover

■ **EASY REMOVER** 27-37



Light, Fast And Efficient

Vakker® Easy Remover is an Automatic crown and bridge remover for the removal of fixed prosthesis and offers genuine comfort for the patient.

Vakker® Easy Remover enables you to successfully carry out crown and bridge removals without stress for the patient.

14

Dental Anesthesia Delivery Device

vakker[®] Dental Practice Helper!

Dental Anesthesia Delivery Device Star Pen & Super Pen

■ Star Pen 27-69

■ Super Pen 27-64

■ Super Pen 27-64-BL

Three adjustable injection speeds

Doctors can choose different speed modes according to habits, making treatment more convenient and faster.

PDL



Low



Hig



15

Dental Microscope

iSee 9000 Dental Microscope

Surgical Microscope

iSee 9000 27-68

Deliver best optical quality and clear surgical field

Integrated 4K High Definition Video System facilitates doctor-patient communication. You can earn patients' trust by sharing visual information with the mand creating comprehensive patient records.

Real-time video or high-definition recorded video for easy team sharing and communication, as well as case sharing and demonstration




Four-in-one integrated 4K module

4K High Definition Video System:

4096 × 2160 ultra-high resolution, 4 times that of 1080P

Present you with with images and videos for crystal-clear details and true-to-life color, which can be saved in the high-speed SD card / U disk.

Wireless foot pedal allows quick and easy control.



Built-in orange filter+ green filter



Orange filter:
Preventing premature curing of composites



Green filter:
Enable you to see tiny nerves and vessels clearly in the surgical environment with blood

16

Teeth Whitening

Prefilled Teeth Whitening Trays

ACTIVE INGREDIENTS

Glycerin, PVP, Water,
Hydrogen Peroxide,
Alcohol, Xylitol,
Sodium polyacrylate,
Ethylcellulose,
Stevioside, Menthol.

PRODUCT BENEFITS

Simple instant whitening
Convenient and Fast
Easy to Use

ITEM FORM

12%HP Trays

NUMBER OF ITEMS

7



Original 5X LED Teeth Whitening Kit



ACTIVE INGREDIENTS

Glycerin, Water (Aqua),
Carbamide Peroxide,
Xylitol, Carbomer,
PEG-6

PRODUCT BENEFITS

10 Min Fast Teeth Whitening
5 LED Light
Intertek certificated for safety

ITEM FORM

18%CP Gel

NUMBER OF ITEMS

3

Ultimate Teeth Whitening Kit

ACTIVE INGREDIENTS

Glycerin, Water (Aqua),
Carbamide Peroxide,
Xylitol, Carbomer,
PEG-6

PRODUCT BENEFITS

Strength Gel
Gels Were by Oratech LLC
28 LED Powerful Whitening Light
The Effectiveness of The Gel
Has Been Certified By Intertek

ITEM FORM

35%CP Gel

NUMBER OF ITEMS

4



Ultimate Teeth Whitening Gel

ACTIVE INGREDIENTS

Glycerin, Water (Aqua),
Carbamide Peroxide,
Xylitol, Carbomer,
PEG-6

PRODUCT BENEFITS

Strength Gel
Gels Were by Oratech LLC
The Effectiveness of The Gel
Has Been Certified By Intertek

ITEM FORM

35%CP Gel

NUMBER OF ITEMS

4

28X LED Light Teeth Whitening Strips Kit

ACTIVE INGREDIENTS

Glycerin, PVP, Water, Hydrogen Peroxide, Alcohol, Xylitol, Sodium polyacrylate, Ethylcellulose, Stevioside, Menthol.

PRODUCT BENEFITS

28 LED Powerful Whitening Light
USB Magnetic Charging
Extra long and wide Whitening Strip
Non-Sensitive Whitening Strips
Long Lasting Whitening
Easy to Carry

ITEM FORM

20%CP Strip

NUMBER OF ITEMS

10



MySmile UV Ultrasonic Retainer Cleaner

DEEP CLEANING

The retainer cleaner uses 45KHz ultrasonic waves to effectively remove food bits, stains, and plaque from all sides. It has 4 UV lamps that sterilize your dental appliances and eliminate odors while cleaning

2 WORKING MODES

The ultrasonic cleaner features 2 cleaning options: I . 5 minutes of ultrasonic cleaning combined with UV light (green indicator light on) and II . 10 minutes of pulsed ultrasonic cleaning with UV light (blue indicator light on). For stubborn stains

VERSATILE CLEANER

With a generous 200ML capacity, the retainer cleaner machine is perfect for all your dental appliances, including retainers, mouthguards, aligners, dentures, and electric toothbrush heads. It also works great for cleaning rings, jewelry, and other accessories



Pro 28X LED Teeth Whitening Kit

ACTIVE INGREDIENTS

Glycerin, Water (Aqua), Carbamide Peroxide, Xylitol, Carbomer, PEG-6

PRODUCT BENEFITS

28 LED Powerful Whitening Light
USB Magnetic Charging
Intertek certificated for safety

ITEM FORM

22%CP Gel

NUMBER OF ITEMS

3



Teeth Whitening Gel Refill Pack

ACTIVE INGREDIENTS

Glycerin, Water (Aqua), Carbamide Peroxide, Xylitol, Carbomer, PEG-6

PRODUCT BENEFITS

10 Min Fast Teeth Whitening
Intertek certificated for safety

ITEM FORM

18%CP Gel

NUMBER OF ITEMS

3



Remineralization Gel

ACTIVE INGREDIENTS

Potassium Nitrate, Carbopol, Deionized Water, Glycerol, PH Adjust Materials, Essence, Edible Pigment.

PRODUCT BENEFITS

Reduce teeth sensitivity fast
Effectively strengthens tooth enamel
Best partner for teeth whitening

ITEM FORM

Gel

NUMBER OF ITEMS

3

5X Teeth Whitening **LED Light**

PRODUCT FEATURES

Cordless

POWER

Battery Operated

WEARTIME

10mins

PRODUCT BENEFITS

5 LED Light
10 Min Fast Teeth Whitening
Removable Teeth Whitening Tray



28X Teeth Whitening **LED Light**

PRODUCT FEATURES

Waterproof & Cordless

POWER

Rechargeable

WEARTIME

10mins

PRODUCT BENEFITS

28 LED Powerful Whitening Light
USB Magnetic Charging

LP221 Cordless Professional Water Flosser

POWER SOURCE

1800mha Battery & Rechargeable

SPECIAL FEATURE

4 Different Modes x 3 Pressure Levels
OLED Display, UVC Light Sanitizer

PRODUCT BENEFITS

Removes Plaque
Fresh Breath

REPLACEMENT TIPS

8

WATER TANK

335ml



LP211 Cordless Water Flosser

POWER SOURCE

2000mha Battery & Rechargeable

SPECIAL FEATURE

Cordless & 5 Different Modes

PRODUCT BENEFITS

Removes Plaque
Deep Cleaning

REPLACEMENT TIPS

8

WATER TANK

260ml



LP233 Pro Electric Oral Irrigator

PRODUCT MODES

NORMAL (1-10 Levels) For daily and deep cleaning
PULSE (1-3 Levels) Massage gums and for people with mild periodontitis
CHILD (1-3 Levels) For sensitive teeth, new users and children under 8 years old

SPECIAL FEATURE

- UVC Sterilizer Light
- Water Heating Optional
- Water Purifying System
- Customizable Pressure

WATER TANK

800ml



DY156 MySmile Electric Toothbrush

AGE RANGE (DESCRIPTION)

Adult

SPECIAL FEATURE

Rechargeable, 5 Modes with 3 Brush Heads

POWER SOURCE

Battery Powered

TOOTHBRUSH HEADS

3 Heads



48,000 VPM Electric Toothbrush Black & White

AGE RANGE (DESCRIPTION)

Adult

SPECIAL FEATURE

Rechargeable, 5 Modes with 6 Brush Heads,
and Travel Case

POWER SOURCE

Battery Powered

TOOTHBRUSH HEADS

6 HEADS

Powerful Retainer Cleaner

ACTIVE INGREDIENTS

Potassium Caraoate, Sodium Carbonate Anhydrous, Sodium Percarbonate, Sodium Bicarbonate, Citric Acid Anhydrous Granular, VP/VA Copolymer, Polyethylene Glycol 180, Sodium Alkyl Sulfate, Sodium Benzoate, Tetra Acetyl Ethylenediamine, Peppermint Oil, FD&C Blue #2 Dye

PRODUCT BENEFITS

Perfect for all dental appliances
Kill 99.99% of odor causing bacteria
Easy to Use

FLAVOR:

Mint

ITEM FORM:

120 Tablets



Powerful Retainer Cleaner

ACTIVE INGREDIENTS

Potassium Caraoate, Sodium Carbonate Anhydrous, Sodium Percarbonate, Sodium Bicarbonate, Citric Acid Anhydrous Granular, VP/VA Copolymer, Polyethylene Glycol 180, Sodium Alkyl Sulfate, Sodium Benzoate, Tetra Acetyl Ethylenediamine, Peppermint Oil, FD&C Blue #2 Dye

PRODUCT BENEFITS

Perfect for all dental appliances
Kill 99.99% of odor causing bacteria
Easy to Use

FLAVOR:

Orange

ITEM FORM:

120 Tablets

Powerful Retainer Cleaner

ACTIVE INGREDIENTS

Potassium Caraoate, Sodium Carbonate Anhydrous, Sodium Percarbonate, Sodium Bicarbonate, Citric Acid Anhydrous Granular, VP/VA Copolymer, Polyethylene Glycol 180, Sodium Alkyl Sulfate, Sodium Benzoate, Tetra Acetyl Ethylenediamine, Peppermint Oil, FD&C Blue #2 Dye

PRODUCT BENEFITS

Perfect for all dental appliances
Kill 99.99% of odor causing bacteria
Easy to Use

FLAVOR:

Mint & Orange

ITEM FORM:

120 Tablets



Travel & Concentrated Mouthwash

Mouthwash

PRODUCT BENEFITS

Kill Bad Breath Germs, Prevents Cavities,
Fluoride Free

Alcohol Free & Portable

100% natural ingredients

FLAVOR

Mint

AGE RANGE

Adult



(for water flosser)



Travel Toothpaste Tablets

ACTIVE INGREDIENTS

Coconut oil, Xylitol, Erythritol, Calcium carbonate, Sodium cocoyl glutamate, Sodium hydrogencarbonate, Guar gum, Menthol, Nano Hydroxyapatite, Nano silver, Magnesium distearate, Cyanidanol, Aloe Barbadosensis Leaf Extract, Stevioside.

PRODUCT BENEFITS

Protection against teeth sensitivity
Prevention of cavities
Enhanced whitening

FLAVOR:

Mint

ITEM FORM:

45 Tablets

MySmile V34 PAP Teeth Whitening Powder

ACTIVE INGREDIENTS

Glycerin, PVP, Water,
Hydrogen Peroxide,
Alcohol, Xylitol,
Sodium polyacrylate,
Ethylcellulose,
Stevioside, Menthol.

FLAVOR:

Mint

AGE RANGE:

Adult

ITEM FORM:

Powder

PRODUCT BENEFITS

Little Foam
Fluoride Free
PAP Whitening



KY-M228 Teeth Whitening Machine

■ KY-M228 27-3-10

Specifications:

Dimension of frame : Height:480mm, Arm length: 900mm, Wide: 1200mm

Weight of machine: 27kg

Input voltage: AC100V~240V, 50/60 HZ

Output voltage: DC18V

Power: 60W

Wavelength of blue light: 400~460nm

Density of blue light: 300~400mw/cm²

Whitening time: 0~15 minutes for one time, two or three times for a treatment



MD887 Teeth Whitening Machine

■ MD887 27-3-5-BL

Specifications:

Broad Spectrum : 450nm~520nm(Blue Light); 600nm~640nm(Red Light);
390nm~410nm(Purple Light)

Voltage: AC:90V~260V

Light Source: 4~5W/pc(Blue Light); 2~3W/pc(Red Light); 1W/pc(Purple Light)

Total Arm Length: 75cm.

Total Height: 120cm.

Net Weight: 5.25kgs.

Total Weight Packing: 7.84kgs.

Box Size: Long:66.5cm; Wide:45.5cm; High:21.5cm.



vakker® means Beautiful in Norwegian, **vakker**® Wishes to Bring Beautiful Smile to ALL Our Customers.

Our mission is to serve the needs of the dental world by providing affordable, quality dental equipment. With customer satisfaction as our top priority, we are dedicated to providing everything you need promptly. All orders ship from our Southern California location. In collaboration with Woodpecker®, we offer a range of innovative products. For any questions, please don't hesitate to contact our team of qualified experts—we're here to help!



8238 Mayten Ave Rancho Cucamonga CA 91730 USA

